


# Hysteroscopy Special



**Purchase a BR Hysteroscope & Receive a Hysteroscopy Sheath at NO COST**



## SUPERIOR VISUALIZATION & DURABILITY:

We can guarantee a brilliant high definition image because we control the construction of your endoscope from start to finish, from optical calculation to glass rod production. BR Endoscopes provide maximum protection against degrading influences of reprocessing.

## FIBER OPTICS & ILLUMINATION:

- Tapered fiber cone
- High transmission wide angle illumination fibers

## LONG LASTING PERFORMANCE:

- Proximal & distal sapphire window
- Medical grade stainless steel tubing
- No coatings
- Laser welded joints
- Hydrophilic gas-filled optical train

## HYSTEROSCOPY SHEATH BENEFITS:

- High Flow Rate
- Compact Slim Line Design
- Smooth Working Channel
- Non-Freezing Stopcocks

---

### Hysteroscopy Scope & Sheath

---

BR970-2930-302 Hysteroscope, autoclavable, Ø 2.9mm, 30°, 302mm

BR975-2730-210 Hysteroscopy Sheath for Ø 2.7/2.9mm 30° scope, 5FR working channel, continuous flow, single piece sheath with replaceable stopcocks, Ø 5.5mm, working length 210mm

---

**PROMO DATES: June 1ST– December 31ST 2017**



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

**BR SURGICAL** • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • [customerservice@brsurgical.com](mailto:customerservice@brsurgical.com)  
Office: 904.642.1366 • Toll Free: 888.642.1366 • Fax: 904.642.1368 • [www.BRSurgical.com](http://www.BRSurgical.com)