

LIGHT SOURCES

LED • Xenon



LED LIGHT SOURCES

BR900-2453

LED 50-Watt Light Source: 4-Port Turret

BR900-2460

LED 100-Watt Light Source: Single Storz port

FEATURES	BR900-2453	BR900-2460
Light Type	50-Watt LED	100-Watt LED
LED Array Life	>50,000 hours	>20,000 hours
Color Temp	6500°K	6000-6500°K
Ports	Turret (Storz, ACMI, Olympus, Wolf)	Single (Storz), other adapters optional
Light Intensity Control	Manual Variable	Auto-Adjust
Dimensions (WxHxD)	7.6" x 5" x 7.75"	14" x 4" x 10"
Weight	6 lbs.	10.8 lbs.



XENON LIGHT SOURCES

BR900-2183

Xenon 180-Watt Light Source: Single Storz port

BR900-2302

Xenon 300-Watt Light Source: 4-Port Turret

FEATURES	BR900-2183	BR900-2302
Lamp	180-Watt Xenon arc	300-Watt Xenon arc
Prorated Lamp Life	650 Hours (average)	
Color Temp	5600°K	
Ports	Single Storz Port	Turret (Storz, ACMI, Olympus, Wolf)
Light Intensity Control	Manual Variable	Automatic Variable
Dimensions (WxHxD)	14" x 4.75" x 13.4"	
Weight	19.8 lbs.	23.2 lbs.



PORTABLE LED LIGHT SOURCE

BR900-7400

PORTABLE LED LIGHT SOURCE

Portable LED light source, includes rechargeable battery and charger, Storz adapter



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com

HEADLIGHTS & ACCESSORIES

Xenon • LED



FIBER OPTIC HEADLIGHTS

BR900-3300 XENON HEADLIGHT
BR900-3310 LED/HALOGEN HEADLIGHT

With bifurcated cable & deluxe headband

- Rotating ring allows for complete control of the spot size without losing intensity
- Sterilizable joystick for position control
- Ratchet crown adjustment
- Foam-backed comfort liners
- Storz adapter

FEATURES	BR900-3300	BR900-3310
Light Source	Xenon	LED or Halogen
Variable Focus Spot	20mm-110mm	20mm-80mm
Included Cable	3.5mm / 9-foot	5mm / 9-foot



COMPACT CART

BR900-9130 Compact Endoscopic Video

Cart, 1 shelf with handle, 4 lockable twin castors Ø 75mm, fixed height 32 ¾"



(*) - Other port adapters available



WIRELESS LED HEADLIGHT

BR900-3191 WIRELESS LED HEADLIGHT

With 2 rechargeable batteries, single bay charger & shoulder carrying bag

- Integrated LED light source and deluxe headband with crown and cranial support
- Lightweight and cable-free
- 4 hours continuous operating time
- Brightness: 30,000 lux
- Color Temperature 4500K
- Color Rendition Index (CRI): 75
- Variable spot 20mm-110mm at 16"
- Multidirectional light adjustment



WIRELESS LED HEADLIGHT

BR900-3190 WIRELESS LED HEADLIGHT

With rechargeable battery and charger

- Integrated LED light source and leather headband
- Lightweight and cable-free
- Brightness: more than 100,000 lux
- 3.5 hours continuous operating time
- Variable spot
- Multidirectional light adjustment
- High intensity pure white light
- Rechargeable Lithium-ion battery included



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com