



# ENT Sinus System

## Sinus Procedure Kits ~ Instruments

*A simple start for your In-Office Balloon Sinus and other Sinus Procedures*

### BR980-9101 Minor Sinus Kit

*Five sterile kits per box, each kit contains the below:*


- 4-ply 17" x 19" Towel/Drape
- 10cc Control Syringe
- Neuro Sponges 1/2" x 3" (10 Patties/Card)
- 2oz Medicine Cups (3 per pack)
- 25ga x 3 1/2" Spinal Needle
- Half Shallow Tray, 24" x 24" Wrap
- Sterilization Pouch



### BR46-90000 Office Balloon Sinus Set

*Each set contains one each of the following products:*

- Vienna Speculum*
- Bayonet Forceps*
- Frontal Seeker*
- Maxillary Seeker*
- Freer Elevator*
- Hartman-Noyes Forceps*
- Frazier Suction Tubes (5FR. & 7FR.)*
- Filter and Paper Sterilization Indicators*
- Sterilization container w/mat 12"x5½"x3"*



### Disposable Supplies a o ai e

BR82-81001 Paper Sterilization Indicators for Mini Container, box of 100

BR82-80003 Paper Filter, rectangle, box of 100, 200mm x 140mm

BR82-14001 Container it Perforated id, flat, nperforated bottom, sil er, 12 x

x 3

OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES


BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com  
Office: 904.642.1366 • Toll Free: 888.642.1366 • Fax: 904.642.1368 • [www.BRSurgical.com](http://www.BRSurgical.com)



# BR46-90001 Myringotomy Mini Set

*Each set contains one each of the following products:*

- *Micro Ear Forceps*
- *Hartmann Ear Specula (2mm, 3mm, 4mm)*
- *Buck Ear Curettes, size 1 (1 sharp / 1 blunt)*
- *Barbara Needle, Curved*
- *Day Ear Hook*
- *Farrell Cotton Applicator*
- *Filter and Paper Sterilization Indicators*
- *Sterilization container w/mat 12"x5½"x3"*



## ENT INSTRUMENTS

<b>BR08-34011</b>	IRIS Scissors, straight, sharp/sharp, round shank, 4½"	<b>BR46-17116</b>	NOYES Nasal Polypus Forceps, serrated, 10x1.2mm, 6½"
<b>BR08-34111</b>	IRIS Scissors, curved, sharp/sharp, round shank, 4½"	<b>BR46-21002</b>	TAKAHASHI Nasal Cutting Forceps, straight, 2.5x10mm bite, 4½"
<b>BR08-37111</b>	STRABISMUS/Baby Metzenbaum Scissors, blunt/blunt, curved, 4¼"	<b>BR46-23782</b>	WEIL-BLAKESLEY Nasal Thru-Cut Forceps, fenestrated jaw, straight, 2.5mm, shaft length 4¾", overall 7½"
<b>BR10-17012</b>	ADSON Dressing Forceps, serrated, extra delicate, 4¾"	<b>BR46-22398</b>	WEIL-BLAKESLEY Pediatric Nasal Forceps, 2.5mm, straight, shaft length 4¾", overall 7½"
<b>BR10-18012</b>	ADSON Tissue Forceps, 1x2 teeth, extra delicate, 4¾"	<b>BR46-29503</b>	FRAZIER Aspiration Cannula, 30°, Ø 3FR., 7½"
<b>BR44-11000</b>	BUCK Ear Curette, straight, blunt, size 0	<b>BR46-29505</b>	FRAZIER Aspiration Cannula, 30°, Ø 5FR., 7½"
<b>BR44-12000</b>	BUCK Ear Curette, angled, blunt, size 0	<b>BR46-29507</b>	FRAZIER Aspiration Cannula, 30°, Ø 7FR., 7½"
<b>BR44-18001</b>	BILLEAU Ear Loop, size 1, 3 mm, 6¼"	<b>BR46-29510</b>	FRAZIER Aspiration Cannula, 30°, Ø 10FR., 7½"
<b>BR44-18002</b>	BILLEAU Ear Loop, size 2, 3.5 mm, 6¼"	<b>BR46-29512</b>	FRAZIER Aspiration Cannula, 30°, Ø 12FR., 7½"
<b>BR44-18010</b>	BILLEAU Micro Ear Loop, size 00, 2 mm, 6¼"	<b>BR46-29903</b>	BARON Suction Tube, Ø 3FR, Ø 1.0 x 75mm
<b>BR44-18011</b>	BILLEAU Micro Ear Loop, size 0, 2.5 mm, 6¼"	<b>BR46-29905</b>	BARON Suction Tube, Ø 5FR, Ø 1.7 x 75mm
<b>BR18-21413</b>	FRAZIER Skin Hook (dura), sharp, 1 prong, 5-1/8"	<b>BR46-29907</b>	BARON Suction Tube, Ø 7FR, Ø 2.0 x 75mm
<b>BR18-21601</b>	JOSEPH Mucous Hooklet, 1 prong, 6¼"	<b>BR46-34300</b>	VON EICKEN (KILLIAN) Suction Tube, 3.0mm, long curve, 5¾"
<b>BR24-104429</b>	WEBSTER Needle Holder, cross serrated jaw, 4¾" TC	<b>BR46-35300</b>	VON EICKEN (KILLIAN) Suction Tube, 3.0mm, short curve, 5¾"
<b>BR24-104490</b>	WEBSTER Needle Holder, delicate, smooth, standard, 5" TC	<b>BR46-41518</b>	FREER Septum Elevator, sharp/blunt, 7"
<b>BR44-35040</b>	Micro Ear Forceps, serrated, 4mm x 0.9mm, 3¼"	<b>BR46-42322</b>	COTTLE Septum Elevator, graduated, 8¾"
<b>BR44-35070</b>	Micro Ear Forceps, serrated, 7mm x 1.2mm, 3¼"	<b>BR46-43418</b>	Maxillary Ostium Seeker, double ended, ball tips 2.2mm and 2.5mm, 7"
<b>BR46-11201</b>	VIENNA Nasal Speculum, small, 5¾"	<b>BR46-43401</b>	Frontal Sinus Seeker, double ended, 8½"
<b>BR46-11202</b>	VIENNA Nasal Speculum, medium, 5¾"	<b>BR44-99610</b>	PHENOL APPLICATOR, (DUBERSTEN Style) stainless, 0.6mm tynes
<b>BR46-16116</b>	HARTMAN Nasal Polypus Forceps, 6"	<b>BR44-20616</b>	LUCAE Bayonet Dressing Forceps, 6¼"
<b>BR44-19002</b>	DAY Ear Hook, small, 2.5mm tip, 6¼"		
<b>BR44-19003</b>	DAY Ear Hook, extra small, 1.5mm tip, 6¼"		

### CPT CODES

<b>31296</b>	Nasal/Sinus endoscopy, surgical, with dilation of frontal sinus ostium (e.g., balloon dilation)	<b>31297</b>	Nasal/Sinus endoscopy, surgical, with dilation of sphenoid sinus ostium (e.g., balloon dilation)
<b>31276</b>	Nasal/Sinus endoscopy, surgical, with frontal sinus exploration, with or without removal of tissue from frontal sinus	<b>31295</b>	Nasal/Sinus endoscopy, surgical, with dilation of maxillary sinus ostium (e.g., balloon dilation), transnasal or via canine fossa
<b>31287</b>	Nasal/Sinus endoscopy, surgical, with sphenoidotomy	<b>31256</b>	Nasal/Sinus endoscopy, surgical, with dilation of maxillary antrostomy
<b>31288</b>	with removal of tissue from the sphenoid sinus	<b>31267</b>	with removal of tissue from maxillary sinus