




SuperCut Iris Scissor

BR Surgical is proud to introduce our new SuperCut/Serrated Iris Scissor hand-made from Germany, featuring a Razor Edge on one side, with Micro Serrations on the other, allowing you to feel the precise cutting action all the way to the tips of the blades .

The feathered “Razor Edge” blade provides the sharpest edge available in any scissor blade, designed to provide a smooth and precise feel, with less traumatic cutting of tissue when used in a surgical procedure.



The “Micro Serrations” are designed to reduce slippage of tissue within the blades of the scissor, and to hold the tissue in the jaws to allow for exact cutting, trimming and excision.

The SuperCut is only recommended for tissue and not for suture or gauze. SuperCut Scissors are designated by Black Ring Handles.

| Item | Description |
|--------------|--|
| BR08-34011SC | IRIS Scissors, round, straight, sharp/sharp, serrated, black handles, 4½" SuperCut |
| BR08-34111SC | IRIS Scissors, round, curved, sharp/sharp, serrated, black handles, 4½" SuperCut |

Call BR Surgical to receive a copy of “Care & Handling of Surgical Instruments”

OVER 20,000 SURGICAL INSTRUMENTS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224

Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.BRSurgical.com





Ceramic Coated SUPERCUT IRIS SCISSOR

BR Surgical introduces their NEW line of **Ceramic Coated SuperCut Iris Scissors** that provide “*State of the Art*” in cutting performance on every procedure. The smooth Ceramic Coating provides less friction on the surface of the scissors, ensuring an effortless “*Glide Action*” of the scissor blades while cutting through tissue. The high performance ultra-hard surface provides better durability against wear and tear on the scissors, extending the life of the instrument due to easier cleaning. Finally, the matte black coating also provides reduced glare and reflection during use.



BR Surgical’s Ultra-Hard Ceramic Coated Surface Benefits:

- 1) **Non Sticking** Cutting Performance
- 2) **Smoother Glide Action** through the tissue
- 3) **Extended Lifetime** of Instrument Performance
- 4) **Reduced Glare and Reflection** during use

| Item | Description |
|--------------|---|
| BR08-34012CC | IRIS Scissors, round, straight, sharp/sharp, serrated, ceramic coated, 4½" SuperCut |
| BR08-34112CC | IRIS Scissors, round, curved, sharp/sharp, serrated, ceramic coated, 4½" SuperCut |

BR Surgical offers Instrument Spray Lube and Instrument Cleaning Brushes to help protect and extend the life of your instruments

Call BR Surgical to receive a copy of “Care & Handling of Surgical Instruments”

OVER 20,000 SURGICAL INSTRUMENTS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224

Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.BRSurgical.com

