



Clean Care Cross Contamination Products


Infection Prevention Products for the ENT Office

Clean Care Cross Contamination Products for the ENT suite help to facilitate a safe and clean working environment for your patients and staff.

CleanAdapt™

An economical solution to provide a new clean connection at the terminal end of the suction tubing between procedures. The CleanAdapt™ system consist of a miniature clean tube with adapter for connecting suction cannulas to the main suction tube and vacuum system for an ENT suction pump and/or treatment cabinet.

CleanAdapt™ is an easy economical way to change the patient end of your suction tubing between procedures.



BR980-9257: CleanAdapt™, Barbed adapter w/ 3" tubing, 25/Box

CleanTract™

A disposable suction pump kit for your ENT cabinet. CleanTract™ limits the potential for cross contamination by providing disposable items for the components in the vacuum system that come in contact or have the potential to come in contact with infectious material.

To insure a high standard in hygiene the components should be physically inspected and changed daily or weekly based on the amount of usage of the system.



BR980-9256: CleanTract™, ENT Disposable Suction Pump Kit, consisting of Air Filter Pump, Suction Tubing 10ft, and Fluid canister 1200CC, 5/box

BR980-9250: AIR Filter for Pump, 1/4" barb x 1/4" barb fittings, 3 per pack

BR980-9243: 1200CC Fluid Canister, 5 per box

BR980-9251: Disposable Suction Tubing, 10 ft., 1/4" ID, 10 per box



*Work with most manufacturers cabinets
*Shown on BR Surgical treatment cabinet

Over 20,000 Surgical Instruments & Procedure Sets • Endoscopic Video Equipment & Related Medical Devices

BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224
Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.bersurgical.com



Clean Care Cross Contamination Products

Infection Prevention Products for the ENT Office

CleanRap™

An easy to use protective barrier designed for an ENT chair. CleanRap™ protects the patient from the risk of cross-contamination from airborne aerosols and body fluids. The covers are designed to fit ergonomically to the chair to eliminate surgeon disruption during an exam or procedure. The snug fit design is achieved by allowing the CleanRap™ to be adjusted to the size of the headrest and back of the chair.

The covers are available in three types so that you can choose the best option for your exam/procedure needs: Headrest Covers, Armrest Covers, and Back Cover.

CleanRap™: simple economical protection for physicians, staff, patients and chairs



BR980-9260: CleanRap™, Armrest cover, 7"x22", 50/Box

BR980-9261: CleanRap™, Headrest cover 16"x28", 25/Box

BR980-9262: CleanRap™, Back/Headrest cover 36"x41", 25/Box

Flexible Endoscope Disposable Sheaths

Sterile and impervious to bacteria and viruses, these sheaths provide the ultimate in patient, staff and scope protection, without distortion.

- Aseptic loading and disposal. Easy on, Easy off – No slip
- Validated barrier to micro organisms as small as 25nms
- Optically clear lens cover
- Latex-free
- Sterile, individually packaged, 20/box
- 5-year shelf life

BR46-15501: Flexible Rhinolaryngoscope Sheath, sterile, disposable, designed for universal fit of Ø 3.4mm to 3.6mm, 20/Box



Over 20,000 Surgical Instruments & Procedure Sets • Endoscopic Video Equipment & Related Medical Devices

BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224
Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.brsurgical.com