



ENT Treatment Cabinet

ENT CABINETS & PROCEDURE ROOM EQUIPMENT



BR900-7506
OTOCART ENT Treatment Cabinet, 4 drawers, built-in pressure/suction pump, with backslash, 30" w x 49³/₄" h x 24¹/₄" d

BR900-7509
OTOSMART ENT Treatment Cabinet, 3 drawers, built-in pressure/suction pump, 35¹/₄" w x 36¹/₂" h x 19³/₄" d



BR900-7510
OPTIMUS ENT Treatment Cabinet, 6 drawers, built-in suction pump, pressure pump and waste bin, 38" w x 38" h x 21¹/₄" d



COMPLETE PROCEDURE ROOMS & SUITES

OTOCART STATION EXAM ROOM

Includes the OtoCart work station and the BR900-75004 Optomic ENT examination chair. The OtoCart Treatment Cabinet is a hard coated aluminum cart with 4 drawers, built-in pressure/suction pump, and a backslash. The ENT exam chair is completely motorized (foot pedal), 330 degree rotation, with Trendelenburg position.

OTOSMART PROCEDURE ROOM

The Otosmart work station is available with a monitor mount and is designed to hold endoscopy equipment, surgical instruments and additional supplies. This procedure room includes the BR900-75007 ENT examination chair which is completely motorized, has programmable memories, chair rotation 380 degrees (190 left / 190 right) and Trendelenburg position.

OPTIMUS PROCEDURE ROOM

The OPTIMUS ENT treatment unit is a fully autonomous mobile work station for the ENT procedure room. The cabinet is designed to allow you to mount the microscope and monitor on the unit, along with space for the endoscopy equipment, surgical instruments and supplies. This provides a small footprint for optimal space utilization. Also in the procedure room is the BR900-75010 ENT examination chair with the following features: completely motorized, programmable memories, electromagnetic brake, synchronizing backrest and Trendelenburg position.