

FLEXIBLE RHINOLARYNGOSCOPE

Superior German Optics in a Well-Balanced & Durable Scope



BR Surgical is pleased to offer four premium-quality Flexible Rhinolaryngoscopes designed for the needs of the physician in daily practice. High resolution and a specially designed multi-lens objective ensures the optimum picture quality you expect from a high-precision scope.

High quality construction and superior optics make this scope a long-lasting value.
Excellent customer service and support make it a smart purchase.

Ergonomic design
for effortless use



Precise motion
for perfect control



Ultra-slim-line
for gentle insertion



High durability
for cost effectiveness



Outstanding optical quality
for brilliant image

| SPECIFICATIONS | Outside Diameter | Working Length | Viewing Direction | Tip Deflection |
|-----------------|------------------|----------------|-------------------|----------------|
| BR946-34260-300 | 3.4 mm | 300 mm | 0° | Up/Down 130° |

CPT CODES

- 31575 Laryngoscopy, flexible fiber optic; diagnostic
- 31576 Laryngoscopy, flexible fiber optic; with biopsy
- 31577 Laryngoscopy, flexible fiber optic; with removal of foreign body
- 31578 Laryngoscopy, flexible fiber optic; with removal of lesions

LIMITED WARRANTY

BR Surgical equipment carries a 12-month warranty from the date of purchase solely on workmanship and all defects of material. Please contact BR Surgical for complete details.



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com

FLEXIBLE RHINOLARYNGOSCOPE

Sinus Scope Video System



FLEXIBLE RHINOLARYNGOSCOPE VIDEO SYSTEM

| | |
|----------------------|--|
| BR900-1101 | Endoscopic Video Camera , 1CCD, including camera head and coupler, with Built-In USB Image Capture System |
| BR900-2460 | LED 100 watt Light Source with single Storz port |
| BR900-4422 | Medical Grade HD Flat Panel Monitor, with DVI and HD, 22" |
| BR900-5050 | Fiber Light Cable , Storz/Storz (Scope/Source) fitting, 3.5mm x 10 ft. |
| BR900-9120 | Mobile Endoscopic Video Cart , with monitor mount and 2 adjustable shelves, 58" x 17 x 16¼" |
| | (optional) LED 50 watt Light source , with 4-port turret (Storz, Olympus, ACMI, Wolf) |
| BR900-2453 | (optional) LED 50 watt Light source , with 4-port turret (Storz, Olympus, ACMI, Wolf) |
| BR900-2150HLS | (optional) Halogen 150 watt Light Source , with Storz port |
| BR900-2183 | (optional) Xenon 180 watt Light Source with Storz port |



FLEXIBLE SCOPES

| | |
|------------------------|--|
| BR946-34260-300 | Flexible Rhinolaryngoscope , Ø 3.4mm, 300mm |
|------------------------|--|

RIGID ENT/SINUS SCOPES

| | |
|-----------------------|--|
| BR946-2700-187 | Sinus Scope , autoclavable, Ø 2.7mm, 0°, 187mm |
| BR946-2730-187 | Sinus Scope , autoclavable, Ø 2.7mm, 30°, 187mm |
| BR946-2770-187 | Sinus Scope , autoclavable, Ø 2.7mm, 70°, 187mm |
| BR946-3000-150 | Sinus Scope , autoclavable, Ø 3mm, 0°, 150mm |
| BR946-3030-150 | Sinus Scope , autoclavable, Ø 3mm, 30°, 150mm |
| BR946-3045-150 | Sinus Scope , autoclavable, Ø 3mm, 45°, 150mm |
| BR946-3070-150 | Sinus Scope , autoclavable, Ø 3mm, 70°, 150mm |
| BR946-4000-175 | Sinus Scope , autoclavable, Ø 4mm, 0°, 175mm |
| BR946-4030-175 | Sinus Scope , autoclavable, Ø 4mm, 30°, 175mm |
| BR946-4045-175 | Sinus Scope , autoclavable, Ø 4mm, 45°, 175mm |
| BR946-4070-175 | Sinus Scope , autoclavable, Ø 4mm, 70°, 175mm |



PORTABLE LIGHT SOURCE

| | |
|-------------------|---|
| BR900-7400 | Portable LED Light Source , Storz adapter, (includes rechargeable battery and charger) |
|-------------------|---|



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com