

# HYSTEROSCOPY SYSTEM

Video System • Scopes • Procedure Kits



## HYSTEROSCOPY VIDEO SYSTEM

<b>BR900-1101</b>	<b>Endoscopic Video Camera</b> , 1CCD, including camera head and coupler, with Built-In USB Image Capture System
<b>BR900-2460</b>	<b>LED 100 watt Light Source</b> with single Storz port
<b>BR900-4422</b>	<b>Medical Grade HD Flat Panel Monitor</b> , with DVI and HD, 22"
<b>BR900-5050</b>	<b>Fiber Light Cable</b> , Storz/Storz (Scope/Source) fitting, 3.5mm x 10 ft.
<b>BR900-9120</b>	<b>Mobile Endoscopic Video Cart</b> , with monitor mount and 2 adjustable shelves, 58" x 17 x 16¾"
<b>BR900-2150HLS</b>	(optional) <b>Halogen 150 watt Light Source</b> , with Storz port
<b>BR900-2183</b>	(optional) <b>Xenon 180 watt Light Source</b> with Storz port



## HYSTEROSCOPES

<b>BR970-2900-302</b>	<b>Hysteroscope</b> , autoclavable, Ø 2.9mm, 0°, 302mm
<b>BR970-2912-302</b>	<b>Hysteroscope</b> , autoclavable, Ø 2.9mm, 12°, 302mm
<b>BR970-2930-302</b>	<b>Hysteroscope</b> , autoclavable, Ø 2.9mm, 30°, 302mm
<b>BR970-4000-302</b>	<b>Hysteroscope</b> , autoclavable, Ø 4mm, 0°, 302mm
<b>BR970-4012-302</b>	<b>Hysteroscope</b> , autoclavable, Ø 4mm, 12°, 302mm
<b>BR970-4030-302</b>	<b>Hysteroscope</b> , autoclavable, Ø 4mm, 30°, 302mm
<b>BR970-4070-302</b>	<b>Hysteroscope</b> , autoclavable, Ø 4mm, 70°, 302mm

## SHEATHS

<b>BR975-2730-210</b>	<b>Hysteroscopy Sheath</b> for Ø 2.7/2.9mm 30° scope, 5FR working channel, continuous flow, single piece sheath with replaceable stopcocks, Ø
<b>BR975-2700-010</b>	<b>Hysteroscopy Diagnostic Sheath</b> for Ø 2.7/2.9mm 0° & 30° scope
<b>BR975-4000-010</b>	<b>Hysteroscopy Diagnostic Sheath</b> for Ø 4mm 0° & 30° scope



High Flow Rate • Compact Slim Line Design • Smooth Working Channel • Non-Freezing Stop-



## FLEXIBLE INSTRUMENTS

	<b>BR980-7010-005 Flexible Scissors</b> , single action jaws, 5FR x 40cm
	<b>BR980-7013-005 Flexible Grasping Forceps</b> , double-action jaws, 5FR x
	<b>BR980-7014-005 Flexible Biopsy Spoon</b> , double-action jaws, 5FR x 40cm

## CPT CODES

<b>52000</b>	Cystoscopy
<b>58563</b>	Endometrial ablation, thermal, with hysteroscopic guidance
<b>58555</b>	Hysteroscopy, diagnostic (separate procedure)
<b>58558</b>	Hysteroscopy, surgical; with sampling (biopsy) of endometrium and/or polypectomy, with or without D&C
<b>58565</b>	Hysteroscopic Sterilization



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • [customerservice@brsurgical.com](mailto:customerservice@brsurgical.com)

# HYSTEROSCOPY Accessories

## Kits & Instruments



### HYSTEROSCOPY KITS • Latex Free



#### BR980-9600

#### Hysteroscopy Sterile Procedure Kit Qty: 5/Box

- |   |                                  |                                |
|---|----------------------------------|--------------------------------|
| 1 BR980-9200 Hysteroscopy Tubing Set            | 1 Needle, 22GA x 3½" Spinal      | 15 Gauze Pads, 4" x 4", 12-ply |
| 1 BR980-9222 Drape (Under Buttock with Pouch)   | 1 Prep, PVP 2-oz. Bottle         | 1 Tray, 2-compartment          |
| 1 BR980-9210 Drape (Video Camera, Elastic Tips) | 1 Bag, Ziplock 6" x 8"           | 1 Cover, Back Table 44" x 76"  |
| 1 Label, white 1" x 3"                          | 2 Syringes, 10cc Control         | 1 Pouch, 14" x 18" Header Bag  |
| 3 O.R. Towels                                   | 4 Swabs, Rayon 8"                |                                |
| 1 Alcohol Prep Pad                              | 1 Gloves, Non-Latex (Medium Sz7) |                                |



#### BR980-9604

#### In-Office Hysteroscopy Pack Qty: 5/Box

- |   |                               |
|---|-------------------------------|
| 1 BR980-9200 Hysteroscopy Tubing Set            | 2 O.R. Towels                 |
| 1 BR980-9222 Drape (Under Buttock with Pouch)   | 1 Alcohol Swabstick           |
| 1 BR980-9210 Drape (Video Camera, Elastic Tips) | 1 Half Shallow Tray           |
| 1 Drape Sheet                                   | 1 Cover, Back Table 44" x 76" |
| 2 PVP 8" Paint Sponge Stick Pouches             | 1 Sterilization Pouch         |
| 2 Leggings                                      |                               |



#### BR980-9601

#### Sterile Draping Set Qty: 5/Box

- |   |                                |
|---|--------------------------------|
| 1 BR980-9200 Hysteroscopy Tubing Set            | 2 White Absorbent Paper Towels |
| 1 BR980-9222 Drape (Under Buttock with Pouch)   | 1 Cover, Back Table 44" x 76"  |
| 1 BR980-9210 Drape (Video Camera, Elastic Tips) | 1 Pouch, 14" x 18" Header Bag  |

BR SURGICAL ALSO OFFERS A FULL LINE OF COATED INSTRUMENTS FOR LEEP PROCEDURES

### ADDITIONAL PROCEDURE KITS

- **BR980-9515: MAJOR Hysteroscopy Procedure Kit**, (5) per box, consisting of hysteroscopy tubing, OB drape w/pouch, endoscopic cleaning brush and 800cc fluid canister
- **BR980-9502: MINOR Hysteroscopy Procedure Kit**, (5) per box, consisting of hysteroscopy tubing, OB drape w/pouch and endoscopic cleaning brush
- **BR980-9520: Hysteroscopy Convenience Set**, (5) per box, consisting of hysteroscopy tubing set, under buttock drape

### ADDITIONAL KIT ACCESSORIES

- **BR980-9200: Hysteroscopy Tubing Set**, sterile, (10) per box
- **BR980-9201: Hysteroscopy Single Line Irrigation Tubing**, sterile, (10) per box, 10 ft.
- **BR980-9222: OB Drape with Pouch**, non-sterile, (5) per box
- **BR980-9210: Camera Drape**, 5" x 96", sterile, (5) per box
- **BR980-9010: Endoscopic Cleaning Brush**, disposable, non-sterile, Ø 1.5mm x 180cm (brush Ø 3.0mm), (10) per package [for cystoscopes & hysteroscopes]
- **BR970-2000-000: Endoscope Sterilization Tray**, with silicone scope holders for rigid scopes, 17½" x 1¼" x 3½"

### STOPCOCK LEVER & CAP

#### BR975-2730-208

Stopcock Lever and Screw Cap replacement for Hysteroscopy Sheath (for BR975-2730-210)



### SHEATH SEALING CAPS

#### BR980-92705

Endoscopic Sheath Sealing Caps, single use, sterile, latex-free, double seal, (10) per pack (for BR975-2730-210, fits 3-10 FR)



### HYSTEROSCOPY RECOMMENDED CHECKLIST

#### GENERAL INSTRUMENTS & NEEDLE EXTENDERS

- **BR16-17025: FOERSTER Forceps**, serrated, straight, 9½"
- **BR16-23126: BOZEMAN Sponge Forceps**, double-curved, 10½"
- **BR70-50425: SCHROEDER-BRAUN Tenaculum Forceps**, straight, 9½"
- **BR70-59640: KEVORKIAN Endocervical Curette**, with basket, 4mm
- **BR70-47008: Needle Extender**, reusable, 3"
- **BR70-47010: Needle Extender**, reusable, 4"
- **BR70-47015: Needle Extender**, reusable, 6"

#### VAGINAL SPECULA

- **BR70-11022: GRAVES Vaginal Speculum**, medium, right side open, 1-3/8" x 4"
- **BR70-11023: GRAVES Vaginal Speculum**, large, right side open, 1½" x 4½"
- **BR70-11024: GRAVES Vaginal Speculum**, extra-large, right side open, 1¾" x 5¼"

#### CERVICAL DILATOR SETS & UTERINE SOUND

- **BR70-43506: HANK Uterine Dilators**, double-ended, set of 6 (sizes 3 - 6.6mm)
- **BR70-43800: PRATT Uterine Dilators**, set of 8, (sizes 13/15 - 41/43FR)
- **BR70-42000: HEGAR Uterine Dilators**, double-ended, set of 8 (sizes 3/4-17/18mm)
- **BR70-58132: SIMS Uterine Sound**, graduated in cm, malleable, silver-plated, 13"



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com