



Fall ENT Endoscope Rebate Offer

INTEGRATED SOLUTION FOR THE ENT OFFICE

RIGID ENDOSCOPES

Take advantage of these special rebates of up to \$350.00 on our high quality endoscopes!

This is a limited time offer, so contact your BR Surgical Rep today!

BR946-4000-175

Sinus Scope, autoclavable,
Ø 4mm, 0°, 175mm

REBATE \$250.00

BR946-4030-175

Sinus Scope, autoclavable,
Ø 4mm, 30°, 175mm

REBATE \$250.00

BR946-3000-150

Sinus Scope, autoclavable,
Ø 3mm, 0°, 150mm

REBATE \$300.00

BR946-3030-150

Sinus Scope, autoclavable,
Ø 3mm, 30°, 150mm

REBATE \$300.00

BR946-2700-140


Sinus Scope, autoclavable,
Ø 2.7mm, 0°, 140mm

REBATE \$350.00

BR946-2730-140

Sinus Scope, autoclavable,
Ø 2.7mm, 30°, 140mm

REBATE \$350.00



SUPERIOR VISUALIZATION & DURABILITY:

We can guarantee a brilliant high definition image because we control the construction of your endoscope from start to finish, from optical calculation to glass rod production. BR Endoscopes provide maximum protection against degrading influences of reprocessing.

FIBER OPTICS & ILLUMINATION:

- Tapered fiber cone = better light amplification
- High transmission wide angle illumination fibers = superior light transmission

LONG LASTING PERFORMANCE:

- Proximal & distal sapphire window = scratch resistant, no glass corrosion
- Medical grade stainless steel tubing = stronger, more durable material
- No coatings = eliminates chance of coatings coming off during reprocessing over time
- Laser welded joints = seals fiberoptic illumination bundles hermetically from the optical train, to avoid humidity and chemical penetration
- Laser welded joints = increase strength between scope body and sheath
- Hydrophilic gas-filled optical train = anti fogging device
- Autoclavable

These rebate offers cannot be combined with any other offer, discount or promotion. Promotion dates October 1, 2021 until December 31, 2021.



Fall ENT Endoscope Rebate Offer

REDEMPTION FORM

GET YOUR REBATE. IT'S EASY!

- Fill out this form...
- Email it along with a copy of your dealer invoice with serial number(s) to customerservice@brsurgical.com
- or mail them to: BR Surgical | Attn: Marketing
3500 Beachwood Court | Suite 107 | Jacksonville, Florida 32224



\$250.00
REBATE

BR946-4000-175
Sinus Scope, autoclavable, Ø 4mm, 0°, 175mm

BR946-4030-175
Sinus Scope, autoclavable, Ø 4mm, 30°, 175mm

\$300.00
REBATE

BR946-3000-150
Sinus Scope, autoclavable, Ø 3mm, 0°, 150mm

BR946-3030-150
Sinus Scope, autoclavable, Ø 3mm, 30°, 150mm

\$350.00
REBATE

BR946-2700-140
Sinus Scope, autoclavable, Ø 2.7mm, 0°, 140mm

BR946-2730-140
Sinus Scope, autoclavable, Ø 2.7mm, 30°, 140mm

Tell us where to send your Rebate Check.

Serial Number(s) _____

Name of Person Ordering _____

Practice Name _____

Address _____

City/State/Zip _____

Phone _____

Dealer/Branch _____

Dealer Representative _____

This offer is valid through 12/31/2021 for our customers inside the USA and does not apply to institutional or government accounts. This promotion only valid with dealer invoice proving purchase. Rebate Redemption Form must be postmarked within 30 days of dealer invoice. Please allow up to 4 weeks for delivery of rebate.