

# ENT/SINUS SYSTEM

## Video System • Scopes



### ENT/SINUS VIDEO SYSTEM

BR900-12003HD	<b>HD Video Camera</b> with Full 1920x1080-HD Video Output, with coupler and with Built-In USB Image Capture and Video Recording System
BR900-2460	<b>LED 100 watt Light Source</b> with single Storz port
BR900-5050	<b>Fiber Light Cable</b> , Storz/Storz (Scope/Source) fitting, 3.5mm x 10 ft.
BR900-4422	<b>Medical Grade HD Flat Panel Monitor</b> , with DVI and HD, 22"
BR900-9120	<b>Mobile Endoscopic Video Cart</b> , with monitor mount and 2 adjustable shelves, 58" x 17 x 16¾"
BR900-2453	(optional) <b>LED 50 watt Light source</b> , with 4-port turret (Storz, Olympus, AC-MACMI, Wolf)
BR900-2183	(optional) <b>Xenon 180 watt Light Source</b> with Storz port
BR900-2302	(optional) <b>Xenon 300 watt Light Source</b> with 4 port turret
BR900-7400	(optional) <b>Portable LED Lightsource</b> , with Storz connector, includes rechargeable battery and charger



### ENT/SINUS SCOPES

BR946-2700-140	<b>Sinus/Arthroscope</b> , autoclavable, Ø 2.7mm, 0º, 140mm
BR946-2730-140	<b>Sinus/Arthroscope</b> , autoclavable, Ø 2.7mm, 30º, 140mm
BR946-2770-140	<b>Sinus/Arthroscope</b> , autoclavable, Ø 2.7mm, 70º, 140mm
BR946-3000-150	<b>Sinus/Arthroscope</b> , autoclavable, Ø 3mm, 0º, 150mm
BR946-3030-150	<b>Sinus/Arthroscope</b> , autoclavable, Ø 3mm, 30º, 150mm
BR946-3045-150	<b>Sinus/Arthroscope</b> , autoclavable, Ø 3mm, 45º, 150mm
BR946-3070-150	<b>Sinus/Arthroscope</b> , autoclavable, Ø 3mm, 70º, 150mm
BR946-4000-175	<b>Sinus/Arthroscope</b> , autoclavable, Ø 4mm, 0º, 175mm
BR946-4030-175	<b>Sinus/Arthroscope</b> , autoclavable, Ø 4mm, 30º, 175mm
BR946-4045-175	<b>Sinus/Arthroscope</b> , autoclavable, Ø 4mm, 45º, 175mm
BR946-4070-175	<b>Sinus/Arthroscope</b> , autoclavable, Ø 4mm, 70º, 175mm
BR946-34260-300	<b>Flexible Rhinolaryngoscope</b> , Ø 3.4mm, 300mm



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • [customerservice@brsurgical.com](mailto:customerservice@brsurgical.com)

# ENT/SINUS SYSTEM

## Instruments



### ENT/SINUS INSTRUMENTS

<b>BR08-34011</b>	<b>IRIS Scissors</b> , straight, sharp/sharp, round shank, 4½"	<b>BR46-17116</b>	<b>NOYES Nasal Polypus Forceps</b> , serrated, 10x1.2mm, 6½"
<b>BR08-34111</b>	<b>IRIS Scissors</b> , curved, sharp/sharp, round shank, 4½"	<b>BR46-21002</b>	<b>TAKAHASHI Nasal Cutting Forceps</b> , straight, 2.5x10mm bite, 4½"
<b>BR08-37111</b>	<b>STRABISMUS/Baby Metzenbaum Scissors</b> , blunt/blunt, curved, 4½"	<b>BR46-23780</b>	<b>WEIL-BLAKESLEY Nasal Thru-Cut Forceps</b> , straight, 2.0mm, shaft length 4¾", overall 7½"
<b>BR08-40516</b>	<b>COTTLE Nasal Scissors</b> , angled 45°, 6¼"	<b>BR46-23792</b>	<b>WEIL-BLAKESLEY Nasal Thru-Cut Forceps</b> , 45° upwards, 2.5mm, shaft length 4¾", overall 7½"
<b>BR10-17012</b>	<b>ADSON Dressing Forceps</b> , serrated, extra delicate, 4¾"	<b>BR46-29503</b>	<b>FRAZIER Aspiration Cannula</b> , 30°, Ø 3FR., 7½"
<b>BR10-17412</b>	<b>ADSON Forceps</b> , serrated, 4¾" TC	<b>BR46-29505</b>	<b>FRAZIER Aspiration Cannula</b> , 30°, Ø 5FR., 7½"
<b>BR10-17512</b>	<b>ADSON Forceps</b> , w/ cross serrated tying platform, 1x2 teeth, 4¾" TC	<b>BR46-29507</b>	<b>FRAZIER Aspiration Cannula</b> , 30°, Ø 7FR., 7½"
<b>BR10-18012</b>	<b>ADSON Tissue Forceps</b> , 1x2 teeth, extra delicate, 4¾"	<b>BR46-29510</b>	<b>FRAZIER Aspiration Cannula</b> , 30°, Ø 10FR., 7½"
<b>BR14-11109</b>	<b>BACKHAUS Towel Forceps</b> , 3½"	<b>BR46-29512</b>	<b>FRAZIER Aspiration Cannula</b> , 30°, Ø 12FR., 7½"
<b>BR14-12111</b>	<b>LORNA Towel Forceps</b> , 4¼"	<b>BR46-29903</b>	<b>BARON Suction Tube</b> , Ø 3FR, Ø 1.0 x 75mm
<b>BR18-21413</b>	<b>FRAZIER Skin Hook</b> (dura), sharp, 1 prong, 5-1/8"	<b>BR46-29905</b>	<b>BARON Suction Tube</b> , Ø 5FR, Ø 1.7 x 75mm
<b>BR18-21418</b>	<b>FRAZIER Skin Hook</b> (dura), sharp, 1 prong, 7"	<b>BR46-29907</b>	<b>BARON Suction Tube</b> , Ø 7FR, Ø 2.0 x 75mm
<b>BR18-21601</b>	<b>JOSEPH Mucous Hooklet</b> , 1 prong, 6¼"	<b>BR46-34300</b>	<b>VON EICKEN (KILLIAN) Suction Tube</b> , 3.0mm, long curve, 5¾"
<b>BR24-104429</b>	<b>WEBSTER Needle Holder</b> , cross serrated jaw, 4¾" TC	<b>BR46-35300</b>	<b>VON EICKEN (KILLIAN) Suction Tube</b> , 3.0mm, short curve, 5¾"
<b>BR24-104490</b>	<b>WEBSTER Needle Holder</b> , delicate, smooth, standard, 5" TC	<b>BR46-41518</b>	<b>FREER Septum Elevator</b> , sharp/blunt, 7"
<b>BR44-20616</b>	<b>LUCAE Bayonet Dressing Forceps</b> , 6¼"	<b>BR46-42322</b>	<b>COTTLE Septum Elevator</b> , graduated, 8¾"
<b>BR44-35040</b>	<b>Micro Ear Forceps</b> , serrated, 4mm x 0.9mm, 3¼"	<b>BR46-43400</b>	<b>Maxillary Sinus Seeker</b> , double ended, 8½"
<b>BR46-11201</b>	<b>VIENNA Nasal Speculum</b> , small, 5¾"	<b>BR46-43401</b>	<b>Frontal Sinus Seeker</b> , double ended, 8½"
<b>BR46-11202</b>	<b>VIENNA Nasal Speculum</b> , medium, 5¾"	<b>BR46-55700</b>	<b>COTTLE Type Skin Hook</b> , angled, 5½"
<b>BR46-12435</b>	<b>COTTLE Nasal Speculum</b> , 35mm, 5¾"		
<b>BR46-16116</b>	<b>HARTMAN Nasal Polypus Forceps</b> , 6"		

### CPT CODES

<b>31296</b>	Nasal/Sinus endoscopy, surgical, with dilation of frontal sinus ostium (e.g., balloon dilation)	<b>31297</b>	Nasal/Sinus endoscopy, surgical, with dilation of sphenoid sinus ostium (e.g., balloon dilation)
<b>31276</b>	Nasal/Sinus endoscopy, surgical, with frontal sinus exploration, with or without removal of tissue from frontal sinus	<b>31295</b>	Nasal/Sinus endoscopy, surgical, with dilation of maxillary sinus ostium (e.g., balloon dilation), transnasal or via canine fossa
<b>31287</b>	Nasal/Sinus endoscopy, surgical, with sphenoidotomy	<b>31256</b>	Nasal/Sinus endoscopy, surgical, with dilation of maxillary antrostomy
<b>31288</b>	with removal of tissue from the sphenoid sinus	<b>31267</b>	with removal of tissue from maxillary sinus



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • [customerservice@brsurgical.com](mailto:customerservice@brsurgical.com)