

FLEXIBLE RHINOLARYNGOSCOPE

Superior German Optics in a Well-Balanced & Durable Scope



BR Surgical is pleased to offer four premium-quality Flexible Rhinolaryngoscopes designed for the needs of the physician in daily practice. High resolution and a specially designed multi-lens objective ensures the optimum picture quality you expect from a high-precision scope.

High quality construction and superior optics make this scope a long-lasting value.
Excellent customer service and support make it a smart purchase.

Ergonomic design
for effortless use



Precise motion
for perfect control



Ultra-slim-line
for gentle insertion



High durability
for cost effectiveness



Outstanding optical quality
for brilliant image

SPECIFICATIONS	Outside Diameter	Working Length	Viewing Direction	Tip Deflection
BR946-34260-300	3.4 mm	300 mm	0°	Up/Down 130°

CPT CODES

- 31575 Laryngoscopy, flexible fiber optic; diagnostic
- 31576 Laryngoscopy, flexible fiber optic; with biopsy
- 31577 Laryngoscopy, flexible fiber optic; with removal of foreign body
- 31578 Laryngoscopy, flexible fiber optic; with removal of lesions

LIMITED WARRANTY

BR Surgical equipment carries a 12-month warranty from the date of purchase solely on workmanship and all defects of material. Please contact BR Surgical for complete details.



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com

FLEXIBLE RHINOLARYNGOSCOPE

Sinus Scope Video System



FLEXIBLE RHINOLARYNGOSCOPE VIDEO SYSTEM

- BR900-12003HD** HD Video Camera with Full 1920x1080-HD Video Output, with coupler and with Built-In USB Image Capture and Video Recording System
- BR900-2460** LED 100 watt Light Source with single Storz port
- BR900-4422** Medical Grade HD Flat Panel Monitor, with DVI and HD, 22"
- BR900-5050** Fiber Light Cable, Storz/Storz (Scope/Source) fitting, 3.5mm x 10 ft.
- BR900-9120** Mobile Endoscopic Video Cart, with monitor mount and 2 adjustable shelves, 58" x 17 x 16¼"
- BR900-2453** (optional) LED 50 watt Light source, with 4-port turret (Storz, Olympus, ACMI, Wolf)
- BR900-2183** (optional) Xenon 180 watt Light Source with Storz port



FLEXIBLE SCOPES

- BR946-34260-300** Flexible Rhinolaryngoscope, Ø 3.4mm, 300mm

RIGID ENT/SINUS SCOPES

- BR946-2700-187** Sinus Scope, autoclavable, Ø 2.7mm, 0°, 187mm
- BR946-2730-187** Sinus Scope, autoclavable, Ø 2.7mm, 30°, 187mm
- BR946-2770-187** Sinus Scope, autoclavable, Ø 2.7mm, 70°, 187mm
- BR946-3000-150** Sinus Scope, autoclavable, Ø 3mm, 0°, 150mm
- BR946-3030-150** Sinus Scope, autoclavable, Ø 3mm, 30°, 150mm
- BR946-3045-150** Sinus Scope, autoclavable, Ø 3mm, 45°, 150mm
- BR946-3070-150** Sinus Scope, autoclavable, Ø 3mm, 70°, 150mm
- BR946-4000-175** Sinus Scope, autoclavable, Ø 4mm, 0°, 175mm
- BR946-4030-175** Sinus Scope, autoclavable, Ø 4mm, 30°, 175mm
- BR946-4045-175** Sinus Scope, autoclavable, Ø 4mm, 45°, 175mm
- BR946-4070-175** Sinus Scope, autoclavable, Ø 4mm, 70°, 175mm



PORTABLE LIGHT SOURCE

- BR900-7400** Portable LED Light Source, Storz adapter, (includes re-chargeable battery and charger)



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com