

HYSTEROSCOPY SYSTEM

Video System • Scopes • Procedure Kits



HYSTEROSCOPY VIDEO SYSTEM

BR900-12003HD	HD Video Camera with Full 1920x1080-HD Video Output, with coupler and with Built-In USB Image Capture and Video Recording System
BR900-2460	LED 100 watt Light Source with single Storz port
BR900-4422	Medical Grade HD Flat Panel Monitor , with DVI and HD, 22"
BR900-5050	Fiber Light Cable , Storz/Storz (Scope/Source) fitting, 3.5mm x 10 ft.
BR900-9120	Mobile Endoscopic Video Cart , with monitor mount and 2 adjustable shelves, 58" x 17 x 16¾"
BR900-2453	(optional) LED 50 watt Light Source w/ 4 port turret (Storz, Olympus, ACMI, Wolf)
BR900-2183	(optional) Xenon 180 watt Light Source with Storz port



HYSTEROSCOPES

BR970-2900-302	Hysteroscope , autoclavable, Ø 2.9mm, 0°, 302mm
BR970-2912-302	Hysteroscope , autoclavable, Ø 2.9mm, 12°, 302mm
BR970-2930-302	Hysteroscope , autoclavable, Ø 2.9mm, 30°, 302mm
BR970-4000-302	Hysteroscope , autoclavable, Ø 4mm, 0°, 302mm
BR970-4012-302	Hysteroscope , autoclavable, Ø 4mm, 12°, 302mm
BR970-4030-302	Hysteroscope , autoclavable, Ø 4mm, 30°, 302mm
BR970-4070-302	Hysteroscope , autoclavable, Ø 4mm, 70°, 302mm

SHEATHS



BR975-2730-210	Hysteroscopy Sheath for Ø 2.7/2.9mm 30° scope, 5FR working channel, continuous flow, single piece sheath with replaceable stopcocks, Ø
BR975-2700-010	Hysteroscopy Diagnostic Sheath for Ø 2.7/2.9mm 0° & 30° scope
BR975-4000-010	Hysteroscopy Diagnostic Sheath for Ø 4mm 0° & 30° scope



High Flow Rate • Compact Slim Line Design • Smooth Working Channel • Non-Freezing Stop-



FLEXIBLE INSTRUMENTS

	BR980-7010-005 Flexible Scissors , single action jaws, 5FR x 40cm
	BR980-7013-005 Flexible Grasping Forceps , double-action jaws, 5FR x
	BR980-7014-005 Flexible Biopsy Spoon , double-action jaws, 5FR x 40cm

CPT CODES

52000	Cystoscopy
58563	Endometrial ablation, thermal, with hysteroscopic guidance
58555	Hysteroscopy, diagnostic (separate procedure)
58558	Hysteroscopy, surgical; with sampling (biopsy) of endometrium and/or polypectomy, with or without D&C
58565	Hysteroscopic Sterilization



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com

HYSTEROSCOPY Accessories

Kits & Instruments



HYSTEROSCOPY KITS • Latex Free



BR980-9600

Hysteroscopy Sterile Procedure Kit Qty: 5/Box

- | | | |
|-------------------------------------------------|----------------------------------|--------------------------------|
| 1 BR980-9200 Hysteroscopy Tubing Set | 1 Needle, 22GA x 3½" Spinal | 15 Gauze Pads, 4" x 4", 12-ply |
| 1 BR980-9222 Drape (Under Buttock with Pouch) | 1 Prep, PVP 2-oz. Bottle | 1 Tray, 2-compartment |
| 1 BR980-9210 Drape (Video Camera, Elastic Tips) | 1 Bag, Ziplock 6" x 8" | 1 Cover, Back Table 44" x 76" |
| 1 Label, white 1" x 3" | 2 Syringes, 10cc Control | 1 Pouch, 14" x 18" Header Bag |
| 3 O.R. Towels | 4 Swabs, Rayon 8" | |
| 1 Alcohol Prep Pad | 1 Gloves, Non-Latex (Medium Sz7) | |



BR980-9604

In-Office Hysteroscopy Pack Qty: 5/Box

- | | |
|-------------------------------------------------|-------------------------------|
| 1 BR980-9200 Hysteroscopy Tubing Set | 2 O.R. Towels |
| 1 BR980-9222 Drape (Under Buttock with Pouch) | 1 Alcohol Swabstick |
| 1 BR980-9210 Drape (Video Camera, Elastic Tips) | 1 Half Shallow Tray |
| 1 Drape Sheet | 1 Cover, Back Table 44" x 76" |
| 2 PVP 8" Paint Sponge Stick Pouches | 1 Sterilization Pouch |
| 2 Leggings | |



BR980-9601

Sterile Draping Set Qty: 5/Box

- | | |
|-------------------------------------------------|--------------------------------|
| 1 BR980-9200 Hysteroscopy Tubing Set | 2 White Absorbent Paper Towels |
| 1 BR980-9222 Drape (Under Buttock with Pouch) | 1 Cover, Back Table 44" x 76" |
| 1 BR980-9210 Drape (Video Camera, Elastic Tips) | 1 Pouch, 14" x 18" Header Bag |

BR SURGICAL ALSO OFFERS A FULL LINE OF COATED INSTRUMENTS FOR LEEP PROCEDURES

ADDITIONAL PROCEDURE KITS

- **BR980-9515: MAJOR Hysteroscopy Procedure Kit**, (5) per box, consisting of hysteroscopy tubing, OB drape w/pouch, endoscopic cleaning brush and 800cc fluid canister
- **BR980-9502: MINOR Hysteroscopy Procedure Kit**, (5) per box, consisting of hysteroscopy tubing, OB drape w/pouch and endoscopic cleaning brush
- **BR980-9520: Hysteroscopy Convenience Set**, (5) per box, consisting of hysteroscopy tubing set, under buttock drape

ADDITIONAL KIT ACCESSORIES

- **BR980-9200: Hysteroscopy Tubing Set**, sterile, (10) per box
- **BR980-9201: Hysteroscopy Single Line Irrigation Tubing**, sterile, (10) per box, 10 ft.
- **BR980-9222: OB Drape with Pouch**, non-sterile, (5) per box
- **BR980-9210: Camera Drape**, 5" x 96", sterile, (5) per box
- **BR980-9010: Endoscopic Cleaning Brush**, disposable, non-sterile, Ø 1.5mm x 180cm (brush Ø 3.0mm), (10) per package [for cystoscopes & hysteroscopes]
- **BR970-2000-000: Endoscope Sterilization Tray**, with silicone scope holders for rigid scopes, 17½" x 1¼" x 3½"

STOPCOCK LEVER & CAP

BR975-2730-208

Stopcock Lever and Screw Cap replacement for Hysteroscopy Sheath (for BR975-2730-210)



SHEATH SEALING CAPS

BR980-92705

Endoscopic Sheath Sealing Caps, single use, sterile, latex-free, double seal, (10) per pack (for BR975-2730-210, fits 3-10 FR)



HYSTEROSCOPY RECOMMENDED CHECKLIST

GENERAL INSTRUMENTS & NEEDLE EXTENDERS

- **BR16-17025: FOERSTER Forceps**, serrated, straight, 9½"
- **BR16-23126: BOZEMAN Sponge Forceps**, double-curved, 10½"
- **BR70-50425: SCHROEDER-BRAUN Tenaculum Forceps**, straight, 9½"
- **BR70-59640: KEVORKIAN Endocervical Curette**, with basket, 4mm
- **BR70-47008: Needle Extender**, reusable, 3"
- **BR70-47010: Needle Extender**, reusable, 4"
- **BR70-47015: Needle Extender**, reusable, 6"

VAGINAL SPECULA

- **BR70-11022: GRAVES Vaginal Speculum**, medium, right side open, 1-3/8" x 4"
- **BR70-11023: GRAVES Vaginal Speculum**, large, right side open, 1½" x 4½"
- **BR70-11024: GRAVES Vaginal Speculum**, extra-large, right side open, 1¾" x 5¼"

CERVICAL DILATOR SETS & UTERINE SOUND

- **BR70-43506: HANK Uterine Dilators**, double-ended, set of 6 (sizes 3 - 6.6mm)
- **BR70-43800: PRATT Uterine Dilators**, set of 8, (sizes 13/15 - 41/43FR)
- **BR70-42000: HEGAR Uterine Dilators**, double-ended, set of 8 (sizes 3/4-17/18mm)
- **BR70-58132: SIMS Uterine Sound**, graduated in cm, malleable, silver-plated, 13"



OVER 20,000 SURGICAL INSTRUMENTS & PROCEDURE SETS • ENDOSCOPIC VIDEO EQUIPMENT & RELATED MEDICAL DEVICES

BR SURGICAL • 3500 Beachwood Court • Suite 107 • Jacksonville, FL 32224 • customerservice@brsurgical.com