




Myringotomy Knife Blade with Handle



BR06-12435

Myringotomy Knife/Blade with Handle,
sterile, spear tip, angled, 5/box

BR06-12436

Myringotomy Knife/Blade with Handle,
sterile, lance tip, angled, 5/box

KEY FEATURES:

- Stainless steel blade provides optimum precision cutting in myringotomy procedures
- Integrated plastic handle provides superior control and no assembly for reusable chuck handle
- 2-directional cutting for versatility
- Lance and Spear cutting blade styles
- Disposable & Sterile

Over 20,000 Surgical Instruments & Procedure Sets • Endoscopic Video Equipment & Related Medical Devices

BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224
Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.brsurgical.com