

**NEW
PRODUCT!**

ELIMINATE CROSS- CONTAMINATION ON YOUR SUCTION TUBES.

Use BR Surgical Frazier
Disposable Suction Tubes.

Avoid Cross – Contamination and Improve Suction Performance

When processing reusable suction tubes there is the potential for blood and organic debris to clog or block the lumen.

Reduce Risk of Infection

Contaminated tubes create the potential for infection in new patients.


Frazier and Baron Design Like the Reusable

Physician doesn't need to change technique during procedures.

Full Breadth of Line

Physician not limited in lumens size or length.

Eliminates cost and difficulty
associated with cleaning and
maintenance of small lumens.



Over 20,000 Surgical Instruments & Procedure Sets • Endoscopic Video Equipment & Related Medical Devices

BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224

Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.brsurgical.com



Disposable Suction Tubes

SINGLE-USE SOLUTION

BR99-29100	BARON-FRAZIER Suction Tube, single-use, White, 30°, Ø 3FR, 80mm, 10/box
BR99-29103	BARON-FRAZIER Suction Tube, single-use, Light Blue, 30°, Ø 5FR, 80mm, 10/box
BR99-29105	BARON-FRAZIER Suction Tube, single-use, Pink, 30°, Ø 7FR, 80mm, 10/box
BR99-29120	FRAZIER Suction Tube, single-use, White, 30°, Ø 3FR, 110mm, 10/box
BR99-29123	FRAZIER Suction Tube, single-use, Light Blue, 30°, Ø 5FR, 110mm, 10/box
BR99-29125	FRAZIER Suction Tube, single-use, Pink, 30°, Ø 7FR, 110mm, 10/box
BR99-29127	FRAZIER Suction Tube, single-use, Lime Green, 30°, Ø 9FR, 110mm, 10/box
BR99-29129	FRAZIER Suction Tube, single-use, Wine, 30°, Ø 11FR, 110mm, 10/box
BR99-29130	FRAZIER Suction Tube, single-use, Green, 30°, Ø 12FR, 110mm, 10/box

Avoid Cross-Contamination



**Organic Debris/Blood in
Reusable Suction Tube**



**New Reusable
Suction Tube**

Over 20,000 Surgical Instruments & Procedure Sets • Endoscopic Video Equipment & Related Medical Devices

BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224
Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.brsurgical.com