




DERMATOLOGY

SPECIALTY CATALOG

SCISSORS • SUPERCUT



SuperCut Scissors are manufactured with a razor edge on one blade and micro serrations on the other blade. The razor edge blade is designed to provide a smooth, precise and less traumatic cutting of tissue and the micro serrated blade is designed to reduce slippage of tissue within the blades of the scissor. The SuperCut is only recommended for tissue not for suture or gauze. **SuperCut Scissors are designated by Black Ring Handles.**

IRIS SCISSORS • SC

BR08-34011SC IRIS Scissor, straight, 4½" SuperCut

BR08-34111SC IRIS Scissor, curved, 4½" SuperCut



BR Surgical's Ultra-Hard Ceramic Coated Surface Benefits:

Non Sticking Cutting Performance • Smoother Glide Action Through the Tissue
Extended Lifetime of Instrument Performance • Reduced Glare and Reflection During Use

IRIS SCISSORS • SC • CERAMIC COATED

BR08-34012CC IRIS Scissor, round, straight, sharp/sharp, serrated, 4½"
SuperCut, ceramic coated

BR08-34112CC IRIS Scissor, round, curved, sharp/sharp, serrated,, 4½"
SuperCut, ceramic coated



METZENBAUM SCISSORS • SC

BR08-28015SC METZENBAUM Scissor, straight, 5¾" SuperCut

BR08-28115SC METZENBAUM Scissor, curved, 5¾" SuperCut



STRABISMUS SCISSORS • SC

BR08-37011SC STRABISMUS/Baby Metzenbaum Scissor, straight, 4¼" SuperCut

BR08-37111SC STRABISMUS/Baby Metzenbaum Scissor, curved, 4¼" SuperCut




STEVENS TENOTOMY SCISSORS • SC

BR08-36011SCS STEVENS Tenotomy Scissor, straight, sharp, 4¼" SuperCut

BR08-36111SCS STEVENS Tenotomy Scissor, curved, sharp, 4¼" SuperCut



SCISSORS • TUNGSTEN CARBIDE



Tungsten Carbide (TC) Scissors are manufactured with Tungsten Carbide Inserts in the blades of the Scissor. Tungsten Carbide is much harder and longer lasting than Stainless Steel and is designed to provide a longer lasting sharp edge on the Scissor blade for many years. **Tungsten Carbide Scissors are designated by Gold Ring Handles.**

IRIS SCISSORS • TC

- BR08-34411 IRIS Scissor, straight, 4½" TC
- BR08-34511 IRIS Scissor, curved, 4½" TC



OR SCISSORS • TC

- BR08-12413 OR Scissor, straight, sharp/sharp, 5" TC
- BR08-12414 OR Scissor, straight, sharp/sharp, 5¾" TC



METZENBAUM SCISSORS • TC

- BR08-28415 METZENBAUM Scissor, straight, blunt, 5¾" TC
- BR08-28515 METZENBAUM Scissor, curved, blunt, 5¾" TC



STRABISMUS SCISSORS • TC

- BR08-37411 STRABISMUS/Baby Metzenbaum Scissor, straight, 4¼" TC
- BR08-37511 STRABISMUS/Baby Metzenbaum Scissor, curved, 4¼" TC



KILNER SCISSOR • TC

- BR08-31412 KILNER Undermining Scissor, curved, non-serrated, flat tips, 4¾" TC



SCISSORS

IRIS SCISSORS

BR08-33810	IRIS Scissor, straight, 4"
BR08-33910	IRIS Scissor, curved, 4"
BR08-34011	IRIS Scissor, straight, 4½"
BR08-34111	IRIS Scissor, curved, 4½"



IRIS RIBBON SCISSORS

BR08-32209	IRIS Scissor, round, straight, 17mm blade, ribbon type, 3¾"
BR08-32309	IRIS Scissor, round, curved, 17mm blade, ribbon type, 3¾"
BR08-32610	IRIS Scissor, round, straight, 22mm blade, ribbon type, 4⅞"
BR08-32710	IRIS Scissor, round, curved, 22mm blade, ribbon type, 4⅞"



OR SCISSORS

BR08-12011	OR Scissor, straight, sharp/sharp, 4½"
BR08-12014	OR Scissor, straight, sharp/sharp, 5½"



METZENBAUM SCISSORS

BR08-28015	METZENBAUM Scissor, straight, blunt, 5¾"
BR08-28115	METZENBAUM Scissor, curved, blunt, 5¾"



STRABISMUS SCISSORS

BR08-37011	STRABISMUS/Baby Metzenbaum Scissor, straight, 4¼"
BR08-37111	STRABISMUS/Baby Metzenbaum Scissor, curved, 4¼"



KILNER SCISSOR

BR08-31012	KILNER Undermining Scissor, curved, non-serrated, flat tips, 4¾"
------------	--



LITTLER DISSECTING SCISSOR

BR08-20512 LITTLER Dissecting Scissor, curved tips with eye for suture, 4¾"



STEVENS TENOTOMY SCISSORS

BR08-36011 STEVENS Tenotomy Scissor, straight, sharp, 4¼"

BR08-36111 STEVENS Tenotomy Scissor, curved, sharp, 4¼"



GRADLE SCISSOR

BR08-34653 GRADLE Eye Suture Scissor, ribbon handle, slightly curved, sharp, 3¾"



SPENCER STITCH SCISSORS

BR08-80209 SPENCER Stitch Scissor, 3½"

BR08-80210 SPENCER Stitch Scissor, 4¼"



STITCH SCISSOR

BR08-84664 Stitch Scissor, angled, 4½"



WILMER CONVERSE SCISSOR

BR08-35010 WILMER CONVERSE Scissor, angled on flat, 4"




NOYES IRIS SCISSORS

BR08-11012 NOYES IRIS Scissor, spring handle, straight, sharp, 4¾"

BR08-11112 NOYES IRIS Scissor, spring handle, curved, sharp, 4¾"




NEEDLE HOLDERS • TUNGSTEN CARBIDE



Tungsten Carbide (TC) Needle Holders are manufactured with Tungsten Carbide Inserts in the jaws of the Needle Holder. Tungsten Carbide is much harder and longer lasting than Stainless Steel and provides sharper teeth which are designed to provide a more precise and secure grip on the needles. Tungsten Carbide Needle Holders are designated by Gold Ring Handles.

WEBSTER NEEDLE HOLDERS • TC

- BR24-104429 WEBSTER Needle Holder, cross serrated jaw, 4¾" TC
BR24-104490 WEBSTER Needle Holder, delicate, smooth, standard, 5" TC



MAYO-HEGAR NEEDLE HOLDERS • TC

- BR24-19414 MAYO-HEGAR Needle Holder, 5" TC
BR24-19416 MAYO-HEGAR Needle Holder, 6" TC




DERF NEEDLE HOLDER • TC

- BR24-11412 DERF Needle Holder, 4¾" TC



CRILE-WOOD NEEDLE HOLDER • TC

- BR24-16415 CRILE-WOOD Needle Holder, 6" TC



OLSEN-HEGAR NEEDLE HOLDERS • TC


- BR24-41412 OLSEN-HEGAR Needle Holder with Suture Scissor, serrated, 4¾" TC
BR24-41414 OLSEN-HEGAR Needle Holder with Suture Scissor, serrated, 5½" TC



NEEDLE HOLDERS

WEBSTER NEEDLE HOLDER


BR24-104480 WEBSTER Needle Holder, smooth, 5"



MAYO-HEGAR NEEDLE HOLDERS


BR24-18014 MAYO-HEGAR Needle Holder, 5"

BR24-18016 MAYO-HEGAR Needle Holder, 6"



DERF NEEDLE HOLDER

BR24-11012 DERF Needle Holder, 4¾"



CRILE-WOOD NEEDLE HOLDER

BR24-16015 CRILE-WOOD Needle Holder, 6"




OLSEN-HEGAR NEEDLE HOLDERS

BR24-41012 OLSEN-HEGAR Needle Holder with Suture Scissor, straight, serrated, 4¾"

BR24-41014 OLSEN-HEGAR Needle Holder with Suture Scissor, straight, serrated, 5½"




TISSUE FORCEPS • TUNGSTEN CARBIDE



Tungsten Carbide (TC) Tissue Forceps are manufactured with **Tungsten Carbide Inserts in the tips of the Forcep**. The TC platform is in front of the 1x2 teeth on the distal tip and is used to provide a secure grip of the needle, suture or tissue. Tungsten Carbide is much harder and longer lasting than Stainless Steel. **Tungsten Carbide Tissue Forceps** are designated by **Gold** on the proximal end of the Forcep.

ADSON FORCEPS


- BR10-17412 ADSON Forcep, serrated, 4¾" TC
BR10-17512 ADSON Forcep, 1x2 teeth, with tying platform, 4¾" TC



TISSUE FORCEPS

ADSON FORCEPS

- BR10-17012 ADSON Dressing Forcep, delicate, serrated, 4¾"
BR10-18012 ADSON Tissue Forcep, delicate, 1x2 teeth, 4¾"
BR10-180941 ADSON TUBIANA Tissue & Suture Forcep, tying platform, 1x2 teeth, 4¾"
BR10-182609 ADSON Forcep, cross serrations, 1x2 teeth, 4¾"
BR10-18613 ADSON Tissue Forcep, fenestrated handle, delicate, 1x2 teeth, 4¾"
BR10-18614 ADSON Tissue Forcep, fenestrated handle, delicate, serrated, 4¾"
BR10-18812 BROWN-ADSON Tissue Forcep, 7x7 side-grasp teeth, 4¾"



BISHOP-HARMON FORCEPS

- BR42-94005 BISHOP HARMON Iris Forcep, straight, serrated, 0.5mm tip, 3½"
BR42-94008 BISHOP HARMON Iris Forcep, straight, serrated, 0.8mm tip, 3½"
BR42-94206 BISHOP HARMON Iris Forcep, straight, 1x2 teeth, 0.3mm tip, 3½"
BR42-94205 BISHOP HARMON Iris Forcep, straight, 1x2 teeth, 0.5mm tip, 3½"
BR42-94208 BISHOP HARMON Iris Forcep, straight, 1x2 teeth, 0.8mm tip, 3½"



CILIA FORCEP

- BR43-42710 Cilia and Suture Forcep, precision fitted diagonal jaws 3mm wide, 3¾"



DESMARRES CHALAZION FORCEPS

- BR43-25409 DESMARRES Chalazion Forcep, small, 3½"
BR43-25609 DESMARRES Chalazion Forcep, medium, 3½"
BR43-25809 DESMARRES Chalazion Forcep, large, 3½"



SEMKEN TISSUE FORCEP

- BR10-27012 SEMKEN Tissue Forcep, 1 x 2 teeth, 5"



SPLINTER FORCEP

- BR10-80209 Splinter Forcep, plain, 3½"



HEMOSTATIC FORCEPS


HARTMAN HEMOSTATIC FORCEPS

- BR12-21010 HARTMAN Hemostatic Forcep, straight, 4"
BR12-21110 HARTMAN Hemostatic Forcep, curved, 4"



HALSTED HEMOSTATIC FORCEPS

- BR12-22012 HALSTED Hemostatic Forcep, straight, 5"
BR12-22112 HALSTED Hemostatic Forcep, curved, 5"



BACKHAUS FORCEPS

- BR14-11109 BACKHAUS Towel Clamp, 3½"
BR14-11113 BACKHAUS Towel Clamp, 5¼"



BLADE HANDLES

BLADE HANDLE

- BR06-10301 Blade Handle No. 3, graduated



SIEGEL BLADE HANDLE

- BR06-10303 SIEGEL Round Blade Handle



BEAVER BLADE HANDLES


- BR06-10910 BEAVER Scalpel Blade Handle, 4"
BR06-10913 BEAVER Scalpel Blade Handle, 5"



RULER

RULER • STAINLESS

- BR02-34015 Ruler, stainless,
w/ mm/inch graduations,
15cm/6"



CURETTES



FOX CURETTES

BR75-12801	FOX Curette (Ø 1mm), oval, 5½"
BR75-12802	FOX Curette (Ø 2mm), oval, 5½"
BR75-12803	FOX Curette (Ø 3mm), oval, 5½"
BR75-12804	FOX Curette (Ø 4mm), oval, 5½"
BR75-12805	FOX Curette (Ø 5mm), oval, 5½"
BR75-12806	FOX Curette (Ø 6mm), oval, 5½"

BR75-12901	FOX Curette (Ø 1mm), round, 5½"
BR75-12902	FOX Curette (Ø 2mm), round, 5½"
BR75-12903	FOX Curette (Ø 3mm), round, 5½"
BR75-12904	FOX Curette (Ø 4mm), round, 5½"
BR75-12905	FOX Curette (Ø 5mm), round, 5½"
BR75-12906	FOX Curette (Ø 6mm), round, 5½"



EXCAVATOR PROBE CURETTE

BR74-50113	Excavator-Probe Curette, double-end, w/hole, 2.0mm - 1.5mm, 5½"
BR74-50115	Excavator-Probe Curette, double-end, w/hole, 2.5mm - 1.5mm, 5½"



EXTRACTORS

SCHAMBERG EXTRACTORS

BR74-11101	SCHAMBERG Comedone Extractor, slightly curved double-ends, square/crimped small loop, 3¾"
BR74-11102	SCHAMBERG Comedone Extractor, slightly curved double-ends, extra fine loops, square/crimped small loop, 3¾"
BR74-11103	SCHAMBERG Comedone Extractor, slightly curved double crimped ends, 3¾"
BR74-11110	SCHAMBERG Comedone Extractor, slightly curved double square ends, 3¾"
BR74-11111	SCHAMBERG Comedone Extractor, double square ends, 3¾"




WALTON COMEDONE EXTRACTOR

BR74-11814	WALTON Comedone Extractor, double-ended, sharp/round, 5½"
------------	---



UNNA COMEDONE EXTRACTORS

BR74-11314	UNNA Comedone Extractor, spoon ends, 5½"
BR74-11514	UNNA Comedone Extractor, round ends, 5½"



SAALFELD COMEDONE EXTRACTOR

BR74-12110	SAALFELD Comedone Extractor, double-ended, sharp/round, 4"
------------	--



ELEVATORS, RETRACTORS & SKIN HOOKS

LOCKE ELEVATOR

BR74-74011	Locke Elevator, small blade, 3x15mm, 4½"
------------	--




DESMARRES RETRACTORS

BR18-23010	DESMARRES Eyelid Retractor, 10mm blade, 6¼"
BR18-23012	DESMARRES Eyelid Retractor, 12mm blade, 6¼"



CROSS ACTION RETRACTORS

BR18-63410	BR Cross Action Retractor, sharp, 4"
BR18-63510	BR Cross Action Retractor, blunt, 4"



COTTLE SKIN HOOK

BR46-55700 COTTLE Type Skin Hook, angled, 5½"



FRAZIER SKIN HOOK

BR18-21413 FRAZIER Dura Hook, sharp, 1 prong, 5⅞"



BR SKIN HOOKS

BR18-20001 BR Skin Hook, single prong, sharp, 6¼"

BR18-20002 BR Skin Hook, double prong, 5mm, sharp, 6¼"



GUTHRIE SKIN HOOKS

BR18-21201 GUTHRIE Double Hook, small, 2.5mm, 4¾"

BR18-21202 GUTHRIE Double Hook, large, 3mm, 4¾"



JOSEPH SKIN HOOKS


BR18-21601 JOSEPH Skin Hook, 1 prong, 6¼"

BR18-21602 JOSEPH Skin Hook, 2 prongs, 2mm, 6¼"

BR18-21605 JOSEPH Skin Hook, 2 prongs, 5mm, 6¼"

BR18-21607 JOSEPH Skin Hook, 2 prongs, 7mm, 6¼"

BR18-21610 JOSEPH Skin Hook, 2 prongs, 10mm, 6¼"



NAIL NIPPERS

NAIL NIPPER

BR74-33014 Nail Nipper, straight, double spring, with lock 5½"



NAIL NIPPERS • BARREL SPRINGS


BR74-33615 Nail Nipper, double action, sharp concave jaw, barrel spring, 6"

BR74-33616 Nail Nipper, double action, sharp straight jaw, barrel spring, 6"



ENGLISH ANVIL NAIL SPLITTER

BR74-33413 English Anvil Pattern Nail Splitter, double spring, 5"



FINGER RING CUTTER

WALTON FINGER RING CUTTER

BR33-17200 WALTON Finger Ring Cutter, 6½"





BR SURGICAL, LLC
3500 Beachwood Court
Suite 107
Jacksonville, Florida 32224
Office: (904) 642-1366
Fax: (904) 642-1368
Toll Free: (888) 642-1366
www.brsurgical.com



OVER 18,000 PATTERNS &
SURGICAL PROCEDURE KITS

brsurgical.com

INSTRUMENT TRAYS

INSTRUMENT TRAY WITH LID

BR83-11113 Instrument Tray with Lid, 8 $\frac{7}{8}$ " x 5" x 2"
BR83-11114 Instrument Tray with Lid, 12 $\frac{1}{2}$ " x 8 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "



BR83-11113



BR83-11114

INSTRUMENT CARE PRODUCTS

BRUSHES

BR82-17004 Instrument Cleaning Brush, stainless steel bristles, autoclaveable, 7 $\frac{1}{4}$ " (Set of 3)
BR82-17005 Instrument Cleaning Brush, nylon bristles, autoclaveable, 7 $\frac{1}{4}$ " (Set of 3)

