

# OP-C16

## HIGH PRECISION AND IMAGE QUALITY



### Professional - Quality imaging

- VESA Mount for Medical Monitor

### Compact Design

- Two-segment arm with ultra-smooth movement
- Allows for work in very tight consultation spaces

### Powerful Illumination

- LED light source
- High light intensity
- Warm light filter

### Supports

- Floor stand
- Wall mount

### Maximum Optical Quality

- Straight binocular with wide field of view

### Effective Magnification

- 5 magnification settings: 1.0x to 34x  
(depending on the combination of objective lens, binocular and eyepiece selected)

### Wide Focus Range

- Objective 250 mm
- Precise adjustment

### Video System

- Compact beam splitter
- Integrated Full HD camera

