

# OP-S9

## Procedures Chair

### COMFORT AND EXCELLENCE FOR DOCTOR AND PATIENT



producing quality

#### **Tilt**

- Visible precision
- 30° tilt memory (vestibular)
- Horizontal position
- Trendelenburg memory

#### **Maximum load capacity**

- 727 pounds
- Robust electric motor system

#### **Electric brake**

- Locks rotation with the press of a button

#### **Exceptional Stability**

- Sturdy construction with a firm and secure backrest
- Stable base with a compact design

#### **Multiple handling options**

- Side-mounted keyboard
- Pedal (optional)
- Handheld remote (optional)

#### **Versatile Headrest**

- Angle and height adjustable
- Quick-release system for easy vestibular testing
- Adjustable for children as well

#### **Comfortable Cushions**

- Range of 5 colors
- Pur Upholstery: plush and comfortable

#### **Great rotation and smooth turning**

- 380° rotation
- Agility for examining both ears
- Precise maneuvers without vibrations

#### **Minimum seat height**

- Minimum 19"
- Maximum 33.5"
- Easy access for patients with reduced mobility or disabilities

#### **Seamless**

- Facilitates cleaning
- Dirt does not penetrate the cushion

