

# OTOSMART PLUS

## Treatment Cabinet

### OPTIMIZE SPACE AND IMPROVE WORK ERGONOMICS



#### New feature

The microscope can be mounted now on the OTOSMART PLUS cabinet

#### Professional - Quality imaging

- VESA Mount for Medical Monitor

#### Easy Disinfection

- Medical Grade Glass, obstacle-free surface for thorough cleaning

#### Integration of Endoscopes

- Integrated supports for rigid and flexible endoscopes, fiber optic cable and camera head

#### Versatility

Space for endoscopy equipment and medical devices

#### Mobility

- 4 omnidirectional wheels  
- Mobility for cleaning and positioning for use

#### Integrated Suction and Pressure

- Automatic On/Off  
- Suction: 70 L/min  
- Pressure: 2 bar

#### Consultation Instruments

- 3 drawers for surgical instruments & supplies

#### Dicards

- Trash Bin  
- Used instrument tray

