



Frazier Disposable Suction Tubes

A SMART SINGLE-USE SOLUTION

Avoid cross-contamination on your suction tubes.



Avoid Cross-Contamination and Improve Suction Performance

When processing reusable suction tubes there is the potential for blood and organic debris to clog or block the lumen.

Reduce Risk of Infection

Contaminated tubes create the potential for infection in new patients.

Frazier Design Like the Reusable

Physician doesn't need to change technique during procedures.

Full Breadth of Line

Not limited in lumens size or length.

Save Time and Money

Eliminates cost and difficulty of cleaning and maintenance of small lumens.

- Color Coded Handles for Easy Identification of French Size
- Single Patient Use
- Ergonomic Handle
- Atraumatic Distal Tip
- Latex Free
- 10 per Box



Over 20,000 Surgical Instruments & Procedure Sets • Endoscopic Video Equipment & Related Medical Devices


BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224

Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.brsurgical.com




Disposable Suction Tubes

SINGLE-USE SOLUTION



FRAZIER SUCTION TUBE | 110mm

BR99-29120	FRAZIER Suction Tube, single-use, White, 30°, Ø 3FR, 110mm, 10/box
BR99-29123	FRAZIER Suction Tube, single-use, Light Blue, 30°, Ø 5FR, 110mm, 10/box
BR99-29125	FRAZIER Suction Tube, single-use, Pink, 30°, Ø 7FR, 110mm, 10/box
BR99-29127	FRAZIER Suction Tube, single-use, Lime Green, 30°, Ø 9FR, 110mm, 10/box
BR99-29129	FRAZIER Suction Tube, single-use, Wine, 30°, Ø 11FR, 110mm, 10/box
BR99-29130	FRAZIER Suction Tube, single-use, Green, 30°, Ø 12FR, 110mm, 10/box



BARON-FRAZIER SUCTION TUBE | 80mm

BR99-29100	BARON-FRAZIER Suction Tube, single-use, White, 30°, Ø 3FR, 80mm, 10/box
BR99-29103	BARON-FRAZIER Suction Tube, single-use, Light Blue, 30°, Ø 5FR, 80mm, 10/box
BR99-29105	BARON-FRAZIER Suction Tube, single-use, Pink, 30°, Ø 7FR, 80mm, 10/box

Over 20,000 Surgical Instruments & Procedure Sets • Endoscopic Video Equipment & Related Medical Devices

BR SURGICAL, LLC • 3500 Beachwood Court • Suite 107 • Jacksonville, Florida 32224
Office: (904) 642-1366 • Fax: (904) 642-1368 • Toll Free: (888) 642-1366 • www.brsurgical.com